


APPLICATION NOTE

Practical Use of STYROS™ HQ Simulated Monolith™ Polymeric In the Separation of OVA

Variants.

The superimposed chromatogram depicted below shows a practical rhythm of scouting the optimal conditions of a separation on a **STYROS™ Simulated Monolith™** column as a general rule.



Chromatogram 1

Commercial OVA on **STYROS™ HQ 4 x 300 mm** and **2.1 x 150 mm** at 1 ml/min (478 and 1,733 cm/hr)

Table 1. Operating parameters.

HPLC System.	Agilent 1100 with thermostatted column compartment and quaternary pump.
Columns	STYROS™ HQ/XH 4 X 300 mm (3.77 ml) and STYROS™ HQ/NB 2.1 X 150 mm (0.52 ml)
Mobile phase.	A: 20 mM Tris, pH=8.2 B: A + 1 M NaCl, pH= 8.2
Flow rates	1 ml/min (1,733 cm/hr and 478 cm/hr of linear velocity)
Gradient	4 to 25 % B in 13.5 cv
Temperature	30°C
Detection	280 nm
Injection volume	4 µl and 15 µl
Pressure Drop	22 bar (320 psi) and 41 bar (595 psi)
Sample:	Commercial OVA (as 98 % pure)

Such sequence allows one to quickly reach the optimal conditions of the separation and then use a longer column for added resolution.

Note that the back pressure of the 300 mm column is similar to a 50 mm Monolith column with substantially higher resolving power.

Such low back pressures make these columns very appropriate for process scale pumps.

The sizes of the columns are not limited either.

Unlike Monoliths with high back pressures the Simulated Monolith™ columns can be made in all sizes including preparative columns with low back pressures.

The Simulated Monolith™ columns are ideal for Simulated Moving Bed Chromatography to generate continuous separation processes.

The high capacity of the resin allows the use of salt at the start of the gradient and an overall lower salt for the complete elution while retaining the shape of the peak. No fronting or tailing is seen.

The benefits that Simulated Monolith™ columns can bring to the process can be summed up as follows:

- Absence of leachables
- High chemical stability
- High physical stability
- Availability in different sizes
- High resolution at low and high flow rates
- Low back pressures
- Tolerant to fast changes of buffer
- High capacity
- Possibility of CIP
- Extended lifetime
- High pressure tolerance
- Availability in most chemistry
- A first step towards process scale separations

A direct comparison with of a similar size Monolith column of 4.6 x 50 mm provides a better picture.

	Simulated Monolith™	Monolith
<i>Maximum operating pressure</i>	3000 psi (21 MPa)	1200 psi (8.2 MPa)
<i>Maximum operating temperature</i>	70 °C	70 °C
<i>Recommended flow rate</i>	Up to 1,740 cm/hr	Up to 540 cm/hr
<i>Maximum flow rate</i>	1,800 cm/hr	720 cm/hr
<i>Typical back pressure at 1ml/min</i>	44 psi	290 psi
<i>Solvent compatibility</i>	All HPLC solvents	Most HPLC solvents
<i>pH range</i>	1-14	2-12
<i>Capacity (BSA)</i>	90 mg/ml	18 mg/ml
<i>Buffer changes</i>	No restrictions	Restricted

